



NOTES :

- 1.0 CUSTOMER
- 2.0 SUBMITTED DRAWING FOR APPROVAL.
- 3.0 STANDARD FEATURES INCLUDED.
- 3.1 CONSTRUCTION, EPOXY COATED MILD STEEL.
- 3.2 FRONT SASH AND SIDES PANELS MADE IN CLEAR TEMPERED GLASS.
- 3.3 POWER SUPPLY 110 V 60 HZ.
- 3.5 ON / OFF SWICHT ON CONTROL PANEL.
- 3.6 STANDARD ACCESORIES.

Designed by Carlos Granados	Checked by	Approved by	Date	Date 3/11/2010	
AIR SCIENCE USA			ECO25 + SINK		
ECO25-Gasemt			Edition	Sheet 1 / 1	