



**NOTES :**

- 1.0 CUSTOMER
- 2.0 SUBMIT DRAWING FOR APPROVAL.
- 3.0 STANDARD FEATURES INCLUDED.
- 3.1 CONSTRUCTION, EPOXY COATED MILD STEEL.
- 3.2 FRONT SASH AND SIDES PANELS MADE IN CLEAR TEMPERED GLASS.
- 3.3 POWER SUPPLY 110 V 60 HZ.
- 3.5 OPTIONAL MONITAIR.
- 4.0 PAINT : YUKX WHITE PPL 51273  
BLACK ONYX COLOR EPOXY RESIN COUNTER UP .

Designed by CFG	Checked by	Approved by	Date 6/1/2011	Date
AIR SCIENCE USA		ECO30 UNIT		
ECO30.GAssem			Edition	Sheet 1 / 1