




REVISION HISTORY				
ZONE	REV	DESCRIPTION	DATE	APPROVED
SIDE PANELS	1	TEMPERED GLASSES	5/17/2017	P CH
METAL HEAD	2	FAN MOD RG190	1/14/2019	P CH
METAL HEAD	3	NEW FAN MOD BI2422063TB2M-8B	6/22/2020	P CH
METAL HEAD	4	INTERNAL LED LIGHTS MOUNTED	6/22/2020	P CH
SIDE PANELS	5	CLEAR ACRYLIC ON SIDE PANELS	8/31/2022	P CH



- NOTES:**
- STANDARD ON ALL UNITS:
  - CONSTRUCTION: EPOXY COATED MILD STEEL SIDES PANELS AND HEAD UNIT.
  - FRONT AND BACK PANELS MADE ON 1/4" CLEAR ACRYLIC.
  - CLEAR ACRYLIC ON SIDE PANELS.
  - CENTRIFUGAL FAN BI2422063TB2M-8B.
  - POWER SUPPLY 110V 60Hz. or 220V 50HZ.
  - MAIN SWITCH ON/OFF.
  - MONITORING, LOW AIRFLOW.
  - INTERNAL LED LIGHTS.
  - W/O



GENERAL TOLERANCES IN MILLIMETERS	
LINEAR & ANGULAR	FORM & POSITION
ISO 2768 V	ISO 2868 L

B-B ( 1 : 10 )

Designed by CFG	Checked by	Approved by	Date 6/17/2014	Date
AIR SCIENCE USA			PURAIR 40XT DUCTLESS FUME HOOD	
Purair 40 XT			Edition 5	Sheet 1 / 1