




REVISION HISTORY				
ZONE	REV	DESCRIPTION	DATE	APPROVED
SIDE PANELS	1	TEMPERED GLASSES	5/18/2017	P CH
METAL HEAD	2	FAN MOD RG190	1/9/2019	P CH
METAL HEAD	3	NEW FAN MOD BI2422063TB2M-8B	6/3/2020	P CH
METAL HEAD	4	INTERNAL LED LIGHTS MOUNTED	6/3/2020	P CH
SIDE PANELS	5	CLEAR ACRYLIC ON SIDE PANELS.	8/15/2022	P CH

NOTES:


- STANDARD ON ALL UNITS:
- 1.0 CONSTRUCTION: EPOXY COATED MILD STEEL SIDES PANELS AND HEAD UNIT.
 - 2.0 FRONT AND BACK PANELS MADE ON 1/4" CLEAR ACRYLIC.
 - 3.0 CLEAR ACRYLIC ON SIDE PANELS.
 - 4.0 CENTRIFUGAL FAN BI2422063TB2M-8B.
 - 5.0 POWER SUPPLY 110V 60Hz. or 220V 50Hz.
 - 6.0 MAIN SWITCH ON/OFF.
 - 7.0 MONITORING, LOW AIRFLOW.
 - 8.0 INTERNAL LED LIGHTS.
 - 9.0 W/O



A-A (1 : 10)



B-B (1 : 10)



BLACK POLYPRO SPILLAGE TRAY 25,00 X 625,00 X 970,00 (OPTIONAL)

GENERAL TOLERANCES IN MILLIMETERS	
LINEAR & ANGULAR	FORM & POSITION
ISO 2768 V	ISO 2768 L

Designed by CFG	Checked by	Approved by	Date 9/21/2012	Date
AIR SCIENCE USA		PURAIR 15XT DUCTLESS FUME HOOD		
		Purair 15 XT	Edition 5	Sheet 1 / 1